


# Covid-19

---




**Espline SARS-CoV-2**



# SARS-CoV-2 structure and assay target



<https://phil.cdc.gov/Details.aspx?pid=23311>



**Assay target =  
NP (Nucleocapsid protein)**

NP is released from membrane and RNA by pretreatment with detergent



NP is measured/detected by sandwich immunoassay

Dinesh et al. (2020). Structural basis of RNA recognition by the SARS-CoV-2 nucleocapsid phosphoprotein bioRxiv <https://dx.doi.org/10.1101/2020.04.02.022194>

# Different types of COVID-19 tests



	Real-time PCR	Antigen test		Antibody test
		CLEIA (Lumipulse SARS-CoV-2 Ag)	Immuno-chromatography assay (Espline SARS-CoV-2)	
Purpose	Detecting the presence of antigen ( = current infection)			Determining past infection (= immunity rates in a community)
Sample type	Nasopharyngeal swab fluid or Saliva	Nasopharyngeal swab fluid or Saliva	Nasopharyngeal swab fluid	Blood (Serum)
Target	Viral RNA	Nucleocapsid protein	Nucleocapsid protein	IgM/IgG antibodies
Sensitivity	★★★	★★★	★★	Depends on kit
Specificity	★★★	★★★	★★★	Depends on kit
Test operation	Instrument-based	Instrument-based	Manual (available for PoC)	Manual (available for PoC)
Sample-to-answer time	~ 2 days	<1hr	<1hr	<1hr
Others		Capable of testing 60 or 120 samples per hour on Fujirebio instrument, Lumipulse G600I or G1200		

## Espline SARS-CoV-2 antigen test kits



- Immuno-chromatography (ICA)
- Sample type: nasopharyngeal swab fluid
- Sample volume: 20uL
- Test result TAT: ~30 minutes
- Recommended in Japanese Guideline




# Product overview



## Simple test procedure




## Easy result interpretation






# Assay principle



# Positioning of Antigen tests

- Surveillance testing versus clinical testing
  - Clinical testing:
    - Symptomatic patients
    - Need for high accuracy and sensitivity
    - Not limited by costs
    - Diagnose delay but you can isolate such a person
  - Surveillance testing
    - Asymptomatic individuals
    - Need for speed
    - Speed of reporting is more important than sensitivity
  
- Disadvantage of sensitive tests
  - Unnecessary quarantine days due to identification in the recovery period who still have the virus but below the infectious threshold
  
- NEED FOR:
  - Need for Speed, frequency of testing, speed of result reporting, and as sensitive as possible (but marginally important)
  - Identifying super spreaders





**Figure 1: Surveillance testing effectiveness depends on frequency.** (A) An example viral load trajectory is shown with LOD thresholds of two tests, and a hypothetical positive test on day 6, two days after peak viral load. 20 other stochastically generated viral loads are shown to highlight trajectory diversity (light grey; see Methods). (B) Relative infectiousness for the viral load shown in panel A pre-test, totaling 31% (blue) and post-isolation, totaling 69% (black). (C) Surveillance programs using tests at LODs of  $10^3$  and  $10^5$  at frequencies indicated were applied to 10,000 individuals trajectories of whom 20% would undergo symptomatic isolation near their peak viral load if they had not been tested and isolated first. Total infectiousness removed during surveillance (colors) and self isolation (hatch) are shown for surveillance as indicated, relative to total infectiousness with no surveillance or self-isolation. (D) The impact of surveillance on the infectiousness of 100 individuals is shown for each surveillance program and no testing, as indicated, with each individual colored by test if their infection was detected during infectiousness (medians, black lines) or colored blue if their infection was missed by surveillance or detected positive *after* their infectious period (medians, blue lines). Units are arbitrary and scaled to the maximum infectiousness of sampled individuals.






## Comparison with qRT-PCR (1/2)

Espline (EL) shows as good a detection rate as qRT-PCR for clinical samples with Ct30 or less.

		qRT-PCR		Total	PPV NPV
		Positive	Negative		
EL	Positive	15	2	17	88%
	Negative	21	68	89	76%
Total		36	70	106	
Concordance		42%	97%		

Drill-down of qRT-PCR positive cases

Espline results for aRT-PCR positive samples




Studied by Toho university  
Under article submission (Ishii et al.)



## Comparison with qRT-PCR (2/2)


Espline shows as good a detection rate as qRT-PCR for clinical samples with Ct30 or less.

Positive rate of Espline for the nasopharyngeal swab samples in UVT (Universal Viral Transport) that were tested positive by qRT-PCR




## ESPLINE SARS-CoV-2 Ag rapid test

- Correlation between ESPLINE SARS-CoV-2 and RT-PCR using 72 Japanese clinical samples:
  - The **concordance rate** between assays was **98%** (44/45 cases) for **negative samples** and **37%** (10/27 cases) for **positive samples**.
  - When **comparison** was **based** on the **number of viral RNA copy numbers** as determined by the one-step RT-PCR method, the **concordance rates** between the two assays was:
    - **83%** (5/6 cases) for **samples with 100 copies or more / test**
    - **50%** (6/12 cases) for **samples with 30 copies or more / test**.



## ESPLINE SARS-CoV-2 Ag rapid test

- Correlation between ESPLINE SARS-CoV-2 and RT-PCR using 124 Japanese administrative samples:
  - The **concordance rate** between assays was **100%** (100/100 cases) for **negative samples** and **67%** (16/24 cases) for **positive samples**. Overall concordance was 94% (116/124)
  - When **comparison** was **based** on the **number of viral RNA copy numbers** as determined by the one-step RT-PCR method, the **concordance rates** between the two assays was:
    - **100%** (12/12 cases) for **samples with 1600 copies or more / test**
    - **93%** (14/15 cases) for **samples with 400 copies or more / test**
    - **83%** (16/18 cases) for **samples with 100 copies or more / test**.
  - Same samples solution was used for both and swab has been submerged in the viral transport media





## Cross reactivity

Espline showed no cross reaction with MERAS-CoV N protein, Influenza virus, or Human corona virus.

Cross reaction study using recombinant N protein of SARS, MARS, HCoV, and Influenza virus

	Antigen	Cross reaction
1	SARS-CoV-2	Yes
2	SARS-CoV	Yes
3	MERAS-CoV	No
4	HCoV-229E	No
5	HCoV-DC43	No
6	HCoV-NL63	No
7	HCoV-HKU1	No
8	Influenza Type A H1N1	No
9	Influenza Type A H3N2	No
10	Influenza Type B	No

Cross reaction study using cultured human corona virus

	Cross reaction	
	Run 1	Run 2
229E (ATCC VR-740)	No	No
OC43 (Amsterdam I)	No	No
OC43 (ATCC-VR1558)	No	No
NL63 (Tokyo/SGH-36)	No	No
HKU1 (Tokyo/SGH-15)	No	No

Conducted by Suzuki T., National Institute of Infectious Diseases

## ESPLINE SARS-CoV-2 Ag rapid test

### Espline SARS-CoV-2



- Immuno-chromatography (ICA)
- Sample type: nasopharyngeal swab fluid
- Sample volume: 20uL
- Test result TAT: ~30 minutes
- Japan PMDA cleared on May 13, 2020
- Reimbursement eligible on May 13, 2020
- Recommended in Japanese Guideline on May 13, and June 16, 2020

- Espline is a qualitative assay

Easy result interpretation



- Espline shows a detection rate as good as qRT-PCR for clinical samples with Ct30 or less.





<https://www.mhlw.go.jp/content/000640554.pdf>

<Unofficial translation by Fujirebio>

### [Guideline for the use of SARS-CoV-2 Antigen detection kit](#)

Initial version: May 13, 2020

Updated: June 16

Ministry of Health, Novel Coronavirus Response Headquarters

This guideline describes the approach and important notes to promote the appropriate use of the rapid antigen detection kit (Product name: Espline SARS-CoV-2, manufacturer: Fujirebio, Inc.) for the diagnosis of Novel Coronavirus infection (COVID-19). At this point the purpose of the use of this kit is to detect positive patients quickly.

When tested positive by this kit, it can be considered as a definitive diagnosis. For a patient with suspected infection and **in between the 2<sup>nd</sup> to 9<sup>th</sup> day after symptom onset, negative result of this test kit does not require additional PCR testing.** On the other hand, this kit requires a higher amount of virus in sample than PCR, at this point **it is not suited to use this kit to a patient who has not presented symptom, because of the detectability limitation of this kit.**



# SD BIOSCIENCE vs ESPLINE

## ■ Sensitivity with recombinant antigen

Ag Concentration (pg/mL)	EL 3×	SD BIO
32,000	4.0	+
8,000	3.0	+
2,000	1.5	+w
1,000	1.0	-
500	1.0	-
250	1w	-
125	-	-
62.5	-	-
31.25	-	-

**Protocol**


Sample: SARS-CoV-2 NP Recombinant Ag

Sample processing :

EL 3×: Sample processing solution 2 10µL + sample 20µL  
 ⇒5min⇒ drop 20µL

SD BIO: Sample processing solution 75µL + sample 15µL  
 ⇒ drop 80µL

- Espline + Sample processing solution 2 showed approx. 10 times higher than SD Bioscience POCT.




## Where is the test being used

---

- Rapid test demand in following areas;
  - Emergency test when PCR is not available (quick results, anywhere)
  - Screening test to identify patients for further test (to save PCR capacity)
  - Test to healthcare professionals
  - Test to broad public; back to work/school/weekly test etc.
- Competitor kits in market
  - Korean (SD Biocensor), Chinese and Japanese (Denka)
  - However, based on internal evaluation, Espline performs x10 higher sensitivity compared to other kits
  - Espline has highest performance for rapid test, not using instrument

# Japanese kit



European kit will not contain a Swab

# Process to CE-marking of ESPLINE SARS-CoV-2

---



# timeline



Due date	Task	responsibility
31 july	Technical file ready (IFU,...)	FRI
11 August	CE-Marking	FRE
11 August	Press release	FRI - FRE
15 September	Kits in warehouse	FRE



# ESPLINE® SARS-CoV-2



P/N	Product name	test/box	Content		Store
260326	ESPLINE® SARS-CoV-2	100 tests	Reaction cassette	10 tests/Box x 10	2~30°C
			Sample Extraction Solution (Squeeze tube)	5 tubes/Bag x 20	
			Applicator Tip	10 tips/Bag x 10	

# FRE support

---

[Victor.ruiz@fujirebio.com](mailto:Victor.ruiz@fujirebio.com)

Marketing Manager

[Veronica.olmedo@fujirebio.com](mailto:Veronica.olmedo@fujirebio.com)

[sales.support@fujirebio.com](mailto:sales.support@fujirebio.com)

[Rikkert.maertens@fujirebio.com](mailto:Rikkert.maertens@fujirebio.com)

Marketing director